

# URINE & SALIVA BASED **DRUG TESTERS**

Alcohol Countermeasure Systems Corp  
60 International Boulevard  
Toronto, Ontario M9W 6J2  
CANADA

[info@acs-corp.com](mailto:info@acs-corp.com)  
[acs-corp.com](http://acs-corp.com)



URINE TESTS

**DISCOVER™  
VISION™**

SALIVA TEST

**REVEAL™**

# URINE & SALIVA BASED DRUG TESTERS



## URINE TESTS



## SALIVA TEST



### DISCOVER™ Urine Drug Test Panel

The DISCOVER single drug panel uses the same accurate strips and testing procedure as the multi-tester, but in a smaller, more convenient size. It is designed for testing one drug at a time.

Available single strips include:

- Amphetamines
- Barbiturates
- Cocaine
- Methamphetamines
- Morphine
- Phencyclidine
- Marijuana

### DISCOVER™ Urine Multi-Drug Tester

DISCOVER multi-drug tester is a convenient urine-based analysis kit which allows detection of up to six drugs at one time. Drug strips can be individually chosen for each test subject and results are viewable in minutes.

| Product | Description Test | Sensitivity |
|---------|------------------|-------------|
| AMP     | Amphetamines     | 1000 ng/ml  |
| BAR     | Barbiturates     | 200 ng/ml   |
| COC     | Cocaine          | 300 ng/ml   |
| METH    | Methamphetamine  | 1000 ng/ml  |
| MOR     | Morphine         | 2000 ng/ml  |
| PCP     | Phencyclidine    | 25 ng/ml    |
| THC     | Marijuana        | 50 ng/ml    |

- FDA 510 (K) cleared for professional use
- SENSITIVE. SAMHSA defined levels
- ACCURATE. 92.9-99.8% specific
- FAST. Results in 4-7 minutes

### VISION™ Urine Multi-Drug Test Cup

The VISION Drugs of Abuse Test Cup is a reliable and accurate method of urine-based drug testing. It's single-step operation allows the operator to secure the lid and initiate the test automatically for up to twelve drugs simultaneously. Results appear within five minutes and remain stable for one hour. A convenient port on the top allows the technician to retrieve the urine sample to conduct further drug analysis without having to remove the lid.

### REVEAL™ Saliva Multi-Drug Tester

| Product | Description Test | Sensitivity |
|---------|------------------|-------------|
| AMP     | Amphetamines     | 50 ng/ml    |
| BZD     | Benzodiazepines  | 20 ng/ml    |
| COC     | Cocaine          | 20 ng/ml    |
| METH    | Methamphetamine  | 50 ng/ml    |
| MTD     | Methadone        | 30 ng/ml    |
| OPI     | Opiates          | 40 ng/ml    |
| PCP     | Phencyclidine    | 10 ng/ml    |
| THC     | Marijuana        | 12 ng/ml    |

- Easy to use
- One step simultaneous testing
- Results in 10 minutes
- Non-invasive saliva testing